




KITTIWAKE
monitoring innovation



ANALEXpql

PQ[®] Technology Ferrous Debris Monitor

The new generation ANALEX_{pqL} ferrous debris monitor adds or upgrades PQ[®] Index measurements for your laboratory.

The ANALEX_{pqL} is a completely new instrument developed to replace the widely used pqM and pqA. Using all new hardware, software and the latest technology the highly accurate ANALEX_{pqL} is an essential addition to the modern used oil analysis laboratory.

The ANALEX_{pqL} ferrous debris monitor provides the most accurate means of detecting and measuring ferrous wear debris in lubricating oils, hydraulic oils and greases.

Benefits

- *Your customers' risks reduced* – Wear debris measured, irrespective of particle size and distribution
- *Your customers' confidence increased* – Established, trusted PQ Index and improved accuracy & repeatability
- *Your laboratory's sample capacity increased* – Measurement time more than halved, sample ID entry simplified and device footprint reduced
- *Your laboratory's uptime increased* - Improved reliability and durability
- *Simple to integrate into your processes* - Sample container flexibility
- *Easy to upgrade to these benefits* - Backwards data compatibility
- *Designed to fit into your laboratory's systems* - PC & Ethernet Connectivity and 4 x USB Ports



Specifications

Repeatability:	±4PQ or ± 1% [whichever is greater]
Resolution:	1 PQ
Cycle Duration:	4 Seconds
Minimum PQ Detection:	5 PQ
Size & Weight:	370mm x 230mm (3.6kg)

Product Code	Description
FG-K17000-KW	ANALEX _{pqL} - PQ Instrument

Optional Extras

FG-K18300-KW	Printer Kit
FG-K18311-KW	Barcode Scanner
FG-K18305-KW	Spare Paper
Robot for sample automation is available upon request	

Consumables

FG-K3-208-KW	100ml HDPE Sample Bottles (288)
FG-K3-207-KW	100ml PET Sample Bottles (288)
FG-K15005-KW	2ml plastic sample pots for oils & greases (3000)



Kittiwake Developments Ltd
 3-6 Thorgate Road Littlehampton
 West Sussex, BN17 7LU United Kingdom
 Tel: +44 (0) 1903 731470
 Email: sales@kittiwake.com
 Web: www.kittiwake.com

Kittiwake GmbH
 Marie-Curie-Str. 5
 25337 Elmshorn Germany
 Tel: +49 (4121) 700 890
 Email: info@kittiwake.de
 Web: www.kittiwake.de

Kittiwake Asia Pacific
 E-8-6 Block E Megan Avenue 1
 189 Jalan Tun Razak Kuala Lumpur 50400
 Malaysia
 Tel: +60 (3) 2333 8906
 Email: zainudiny@kittiwake.com
 Web: www.kittiwake.com



Kittiwake Proactive Technologies Pvt Ltd
 405 Ansals Majestic Tower 17
 Block G1 Vikas Puri Community Centre
 New Delhi 11008 India
 Tel: (+91) 11 65356292
 Email: deepaks@kittiwake.com

Kittiwake Americas
 5177 Richmond Avenue Suite 1145
 Houston TX 77056
 Tel: +1 713 255 7255
 Email: keithm@kittiwake.com
 www.kittiwake-americas.com